# BEAUTIFLY

Lumix  
P R O  
series

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- [www.beautify.eu](http://www.beautify.eu)

- [www.youtube.com/beautifyeu](https://www.youtube.com/beautifyeu)

Size	56x187x45.8 mm
Weight	about 154 g
Rated power	about 10W
Rated voltage	DC 5V
Rated current	about 2A
Timer	6 min.
Working mode	RF (1000 KHz) / EMS 5-9 $\mu$ A
Battery	7.4V 600mAh

**Safety:** Before using the product, please read this user manual and follow the instructions contained in it during use.

Device description [Page 2-f1]:

1. Power button / 2. RF function indicator / 3. EMS function indicator / 4. Charging base

User manual:

**NOTE!** Before using the device for the first time, it is recommended to perform a skin sensitivity test. The test should be performed 24 hours before using the Lumix PRO for the first time. To do this, apply a small amount of gel/serum (or other cosmetic with a slippery consistency) to the inner part of the arm and perform the treatment on this area. Make sure that the skin does not react negatively to the gel/serum (the cosmetic you used) or the treatment. It is recommended to set the power to low for the first treatment session and then adjust it for subsequent treatments.

1. Wash your face.
2. Apply an appropriate amount of gel/serum or other cosmetic with a slippery consistency to the area to be used (face or neck).
3. Long press the device button to turn it on.
4. Select the preferred mode.
5. Divide the face into 4 treatment areas: right cheek, left cheek, forehead, outer corners of the eyes.
6. Place the device on the skin, remembering to press the device lightly to ensure full contact of the head with the skin.
7. Start moving the device in slow, very small circular movements on a small part of the treatment area (about twice the size of the treatment head). The device's energy will be released into the skin, providing comfort and warmth to the skin.

**NOTE:** Do not hold the device in one position without moving for too long.

8. Continue working over the entire treatment area in slow circular motions for 2-4 minutes. Make sure to maintain uniform energy by massaging the entire treatment area.

**NOTE:** The treatment effect can be enhanced by performing "lifting" movements from the jawline upwards towards the cheeks.

9. After completing the 2-4 minutes of treatment, you can move on to the next treatment area by repeating steps 6-8. In RF + EMS mode (see the "Modes" section of the instructions), it is recommended to start with a low power level.

**Important:** use RF + EMS mode **ONLY** on the jawline area. Do not use RF + EMS mode on other areas of the face.

10. Begin moving the device in very small linear movements along the jawline.

**NOTE:** For optimal results, use the device for 4 minutes on each side. After massaging the entire treatment area, the treatment session is finished. Natural effects of the treatment may include skin redness (erythema) and slight swelling (puffiness) in the treatment area. Both should disappear after a short time.

11. Long press the power button to turn off the device.

12. Clean the energy head of the device with a damp cloth or paper towel, then place the device on the protective base to protect the energy head until the next use.

Treatment recommendations:

- The treatment of the entire face takes 8-12 minutes. To shorten the treatment time, focus only on the selected areas you want to treat. Treatment times vary for each person.
- For optimal results, use the Lumix PRO 2-3 times a week for 6-8 weeks.
- Allow at least 1 day between treatments on the same treatment area.
- To maintain long-lasting effects, after completing the full course, we recommend using Lumix PRO 2-4 times a month or as needed.

Modes:

In start-up mode, short press the button to switch between Mode 1 - Mode 5. The device will make a

sound and the LED indicator color will change:

Mode 1: RF treatment (level 1), the LED indicator will be steady green.

Mode 2: RF treatment (level 2), the LED indicator will be steady blue.

Mode 3: EMS treatment (level 1), the LED indicator will be steady green.

Mode 4: EMS treatment (level 2), the LED indicator will be steady blue.

Mode 5: RF and EMS treatment (the LED indicator will be steady blue when automatically switching between RF and EMS therapy).

Temperature sensor: The treatment heads are equipped with a precise NTC temperature sensor.

After the RF treatment starts, the sensor will detect the skin temperature, if the temperature is lower than 42°C, the RF treatment will start working. If it is above 42°C, the RF treatment will be stopped.

Automatic timing: The automatic timing function sets the time to 6 minutes. After the time has elapsed, the product stops working and must be restarted to continue using it.

Charging:

1. For the first time use, if the device is not fully charged, please connect it to the charger.
  2. Place the device on the charging base, then insert the USB connector into the charger plug and connect it to the power supply. The LED indicator will be red during charging, and the red light will go out after full charge.
  3. During charging, the 415nm wavelength lamp will disinfect the treatment heads for 5 minutes.
  4. The charging time is about 120 minutes, and the working time of the device is about 50 minutes.
- NOTE: When the device's battery is too low, the buzzer will sound every 5 seconds.

Warnings and recommendations:

1. The device is not waterproof, so do not wash it or place it near water, and keep it in a dry place out of the reach of children.
2. Before each use, check that your hands are dry and the device and power adapter are not damaged; use only the supplied cable and power adapter.
3. Always turn off the device after use or before cleaning it.
4. Use a gel or serum during the treatment to ensure proper glide and effectiveness of the RF/EMS, as well as better skin hydration.
5. The device has an automatic temperature regulation mechanism that protects the skin from overheating. Use the device only according to the instructions, avoiding use on the eyelids, around the eyes, thyroid, lips, ears, breasts and genitals.
6. Use the device comfortably, and in the event of any discomfort, stop the treatment immediately.
7. The device is not a toy. Do not use the device if you are under 18 years of age, have a pacemaker, internal defibrillator, metal implants in the treatment area (except dental fillings), currently have or have had cancer (especially skin cancer) or have pre-cancerous moles. Do not use the device if you have a serious medical condition such as heart problems, are pregnant or breastfeeding, have a weakened immune system due to illness (e.g. AIDS) or immunosuppressive drugs, or have heat-induced diseases such as recurrent herpes. Avoid using the device if you are prone to keloid scars, have difficulty healing wounds or have delicate skin, or if you have had surgery, laser or peeling in the treatment area in the last three months and are still healing. If you have permanent fillers in the treatment area or temporary fillers, consult your doctor before using the device, as well as if you have had Botox in the last three months. Do not use the device on areas with suspected skin lesions or if you are taking medications, herbs, vitamins or supplements that may increase skin sensitivity or dryness.

More information about the company and products [page 2 - f2]:

- [www.beautify.eu](http://www.beautify.eu)

- [www.youtube.com/beautifyeu](https://www.youtube.com/beautifyeu)

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**[EN] Rules for handling electrical and electronic equipment:** Used electrical and electronic equipment and batteries and accumulators must not be mixed with other waste. This is indicated by the crossed-out bin symbol on the product packaging. Electrical and electronic equipment and batteries and accumulators may contain hazardous substances, mixtures and components that are harmful to the environment and human health. Improper handling or damage to them may result in damage to health or environmental pollution.

Proper handling of used electrical and electronic equipment and batteries and accumulators in the household contributes positively to the reuse and recovery, including recycling, of used equipment.

Used electrical and electronic equipment and batteries and accumulators can be returned to local selective waste collection points. The distributor is obliged to collect used equipment from households free of charge at the point of sale, provided that the used equipment is of the same type and has performed the same functions as the equipment sold. A distributor, when delivering household equipment to a purchaser, is obliged to take back used household equipment free of charge at the place of delivery of such equipment, provided that the used equipment is of the same type and has performed the same functions as the equipment delivered. A distributor running a retail outlet with a sales area within the meaning of art. 2 item 19 of the Act of 27 March 2003 on spatial planning and development of at least 400 m<sup>2</sup> devoted to the sale of household equipment, is obliged to accept, free of charge, in that outlet or in its immediate vicinity, used household equipment, none of the external dimensions of which exceeds 25 cm, without the need to purchase new household equipment.

The product was introduced to the market after 13 August 2005.

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**[EN] Summary declaration of conformity:** The manufacturer hereby declares under its sole responsibility that the above product complies with Directive 2014/30/EU of the European Parliament and of the Council of 26 February 2014 on the harmonisation of the laws of the Member States relating to electromagnetic compatibility and Directive 2011/65/EU of the European Parliament and of the Council of 8 June 2011 on the restriction of the use of certain hazardous substances in electrical and electronic equipment, together with Commission Delegated Directive (EU) 2015/863 of 31 March 2015 amending Annex II to Directive 2011/65/EU of the European Parliament and of the Council as regards the list of restricted substances.

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