PL EN DE FR ES PT IT CZ SK SL EE LT LV

NB FI HU RO AR RU EL HR SV DA NL TR BG

# BEAUTIFLY

Skin

# Body

series

## Slimming anti-cellulite massager

/beautifly.eu

@beautifly.eu

beautifly.eu

/beautiflyeu

@beautifly.eu

f1

f2

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**Technical data:**

Rated power:	16W
Rated voltage:	9V
Material:	ABS
Net weight:	~0.8kg

**Security**

Before using the product, read this manual and follow the recommendations contained therein when using the device.

**Cellulite and Body Shaping**

The Skinn device is designed to effectively fight cellulite of any type and for body shaping. Regular massage sessions using it break down fat cells, improve blood and lymph circulation, eliminating cellulite and preventing its formation. The Skinn massager makes the skin smooth and firm.

**Application of Red Light Therapy**

Red light therapy available in the device helps to firm the skin, improve its elasticity, reduce fat tissue and relieve muscle tension. It is also helpful in reducing spasms and swelling, offering pain-relieving and relaxing effects.

**Safety rules**

1. Charging of the device is indicated by the indicator light changing from red to green. After charging, disconnect the charger from the device and the network.
2. Only use the device after disconnecting the charging plug.
3. Avoid using the massager in damp places. Do not immerse it in water or other liquids.
4. Before cleaning the device, disconnect the power supply.
5. Make sure the plug is properly connected to prevent device failure.
6. When connecting the power plug, do not pull the cord.
7. Avoid frequently plugging and unplugging the power cord.
8. If the device is damaged, have it repaired by an authorized service center.
9. The device should not be used by children and people with limited ability to use it independently.
10. The storage place of the device should be inaccessible to children and disabled people.

**Warnings**

People with pacemakers, other metal implants, heart or spine diseases, after surgery, who are pregnant or have recently given birth should not use the device.

**Massage Instructions**

- Arms: Massage both sides evenly from elbows to armpits in circular motions.
- Waist and Belly: Massage the waist up and down and the stomach in a clockwise direction around the navel.
- Legs: Start at the ankles, massaging upwards in straight and then circular motions towards the heart.
- Hips and Buttocks: Massage with movements from the bottom up to reduce cellulite and firm the skin.
- Other Body Parts: Perform the massage in circular or straight movements, always aiming at the heart.

Overlay description [page 2 - f1]:

1. DIVIDED OVERLAY WITH ROUND NUTS - prevents sagging, stimulates blood flow, relieves muscle tension. It resembles lymphatic drainage and its action is mainly to stimulate blood and lymph circulation.
2. BALL PAD - stimulates and firms the muscles and skin, giving it a smoother appearance. It improves circulation in the abdominal area, eliminates toxins and prevents constipation. It has anti-cellulite properties and breaks down fat lumps.
3. FLAT TOPPER – perfect for a relaxing massage with oils
4. ROSETTE OVERLAY – used on thighs and buttocks, improves their elasticity and blood flow.
5. TRIPLE WAVE OVERLAY – relaxes and relieves muscle fatigue
6. MICRO FIT OVERLAY – stimulates collagen, improves hydration, improves skin density and elasticity. Thanks to stimulation, the creams used penetrate deeply into the skin.
7. WAVED OVERLAY WITH LONGITUDINAL AND ROUND PUBS - eliminates cellulite by breaking down fat cells and removing excess water from the body.
8. WAVED OVERLAY – simulates a professional massage toning and improving blood flow
9. MESH COVER – serves to cover the massager cover and protects against fabrics getting entangled in the mechanism.

More information about the company and products [page 2 - f2]

- [www.beautifly.eu](http://www.beautifly.eu)
- [www.youtube.com/beautifly.eu](http://www.youtube.com/beautifly.eu)

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**[EN] Rules for handling electrical and electronic equipment:** Electrical and electronic equipment as well as batteries and accumulators must not be combined with other household waste. This is indicated by the crossed-out wheeled bin symbol on the packaging. Dispose of waste electrical and electronic equipment as well as batteries and accumulators at your local separate collection point or to your dealer. Electrical and electronic equipment and batteries and accumulators may contain hazardous substances, mixtures and components that are harmful to the environment and human health. Incorrect handling or damage to them may lead to damage to health or environmental pollution when the equipment is further disposed of or recycled. The device was placed on the market after August 13, 2005

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**[EN] Summary declaration of conformity:** The manufacturer hereby declares under its sole responsibility that the above product complies with Directive 2014/30/EU of the European Parliament and of the Council of 26 February 2014 on the harmonisation of the laws of the Member States relating to electromagnetic compatibility and Directive 2011/65/EU of the European Parliament and of the Council of 8 June 2011 on the restriction of the use of certain hazardous substances in electrical and electronic equipment, together with Commission Delegated Directive (EU) 2015/863 of 31 March 2015 amending Annex II to Directive 2011/65/EU of the European Parliament and of the Council as regards the list of restricted substances.

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## Shortened warranty terms

This Guarantee is granted by a Polish company, Beautifly Sp. z o. o. based in Warsaw. Plac Bankowy 2, 00-095 Warsaw, Poland. You can find us in the official Polish register of enterprises under the KRS number: 0000887306, NIP: 5252854218 and REGON: 388382021 We guarantee that the equipment you have purchased from us (hereinafter referred to as the "Product") will work properly. If something goes wrong for a reason inherent in the Product, this warranty protects you. The warranty starts from the day you receive the Product and lasts:

24 months, if you are a consumer or a natural person running a business and the purchase of the Product is not of a professional nature for you. Otherwise, the warranty period is 12 months. If the law in your area requires a longer warranty, the warranty period will be extended to meet local law requirements.

In addition, if the Product does not comply with the contract, you have the right to legal protection from and at the expense of the seller. This warranty does not affect these rights - it does not limit, suspend or exclude them.

Detailed warranty conditions, how to submit complaints, what happens when you find a defect, what are exceptions to the warranty, when the warranty expires or is extended, and what are the rules regarding the processing of personal data, you can find on our website: [www.beautifly.eu/warranty-beautifly](http://www.beautifly.eu/warranty-beautifly)