BEAUTIFLY

Longa
PRO
series

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Security:

Before using the product, read this user manual and follow the recommendations contained therein when using the device.

Device description [page 2-f1]:

1. Power
2. RF
3. Laser

The device uses radio wave and laser technology to stimulate hair growth.

Main features:

- 10 laser heads to stimulate hair follicles.
- RF technology to stimulate the scalp.
- LLLT Technology
- Operation via vibration and voice command.
- Anti-scratch design.
- Intelligent timing function allows you to care for your hair in just 15 minutes.

Use:

- Hair loss
- Alopecia areata.
- Androgenetic alopecia.
- Hair weakened due to e.g. disease, anesthesia, pregnancy.

How to use:

1. The Longa hair growth laser should be used in accordance with the action program, i.e. 15 minutes, 3 times a week, preferably every other day, e.g. Monday, Wednesday, Friday. It is recommended to use for another 6 months.
2. The scalp should be freshly washed and the hair should be dry.
3. Long hair should be used before applying styling products because they may obstruct the flow of laser light.
4. Press and hold the power button to turn on your device. Then press the RF radio wave and Laser function to start the treatment. Place the Longa at the roots of your hair and move it from the top

of your head to the nape of your neck (like combing your hair). Repeat the action for 15 minutes. The automatic timer will turn off the device itself after the treatment is completed.

5. After the treatment, clean the head with a soft, slightly damp cotton cloth. Do not use strong cleaning detergents on any part of the device.

Device operation:

1. Power on: Press and hold the power button to turn on, the power button will light up.
2. RF function: in the power-on state, press the RF button to turn on the RF when the RF indicator is on; press it again to turn off RF.
3. Laser function: in the on state, press the laser button to turn on the laser, when the laser pointer is on; press it again to turn off the laser function.
4. Power off: After powering on, press and hold the power button to turn off the laser comb and the power indicator will turn off.
5. Charging: Connect the included USB cable to the power adapter to charge your device. While charging, the power button flashes red. When fully charged, it glows steadily. (For safety reasons, all device functions do not work while charging.)

Securing:

Cleaning and maintenance

1. When the device is dirty, use a dry soft cloth or towel to wipe it.
2. After the treatment, clean the head with a soft, slightly damp cotton cloth. Do not use strong cleaning detergents on any part of the device.
3. Do not immerse the device in water or other liquids. The device is not waterproof.

Storage

1. Keep away from children.
2. Do not place it in a place exposed to direct sunlight, high temperature and humidity.
3. Store in a dry and well-ventilated place.
4. Do not disassemble, repair or modify the device.

Attention:

1. Do not look directly at the red light or the red light in the mirror, which may temporarily irritate your eyes. In addition, do not wear corrective glasses when using the Longa device, otherwise it may cause temporary eye irritation.
2. Do not put Longa into water to prevent electric shock. The warranty does not cover mechanical damage.
3. Do not use Longa device near heat source accessories, which may lead to failure of normal operation.
4. Laser - Avoid direct contact with eyes.
5. If you are taking any medications or products that may cause you to be concerned about light, consult your doctor first.

6. Caution should be exercised in women who are breastfeeding or menstruating.
7. Do not disassemble, modify or transform the product without permission.
8. Do not hit the product as it may cause internal damage.
9. Pregnant women, people with sensitive skin, heart disease, abnormal blood pressure, cancer, patients with cerebrovascular disease, acute illness or other people receiving medical treatment are not allowed to use this product.
10. Product disposal is subject to applicable environmental protection regulations. Otherwise, environmental pollution may occur.

More information about the company and products [page 2 - f2]:

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[EN] Rules for handling electrical and electronic equipment: Electrical and electronic equipment as well as batteries and accumulators must not be combined with other household waste. This is indicated by the crossed-out wheeled bin symbol on the packaging. Dispose of waste electrical and electronic equipment as well as batteries and accumulators at your local separate collection point or to your dealer. Electrical and electronic equipment and batteries and accumulators may contain hazardous substances, mixtures and components that are harmful to the environment and human health. Incorrect handling or damage to them may lead to damage to health or environmental pollution when the equipment is further disposed of or recycled. The device was placed on the market after August 13, 2005

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[EN] Summary declaration of conformity: The manufacturer hereby declares under its sole responsibility that the above product complies with Directive 2014/30/EU of the European Parliament and of the Council of 26 February 2014 on the harmonisation of the laws of the Member States relating to electromagnetic compatibility and Directive 2011/65/EU of the European Parliament and of the Council of 8 June 2011 on the restriction of the use of certain hazardous substances in electrical and electronic equipment, together with Commission Delegated Directive (EU) 2015/863 of 31 March 2015 amending Annex II to Directive 2011/65/EU of the European Parliament and of the Council as regards the list of restricted substances.

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