# BEAUTIFLY

B-Shine

# Blush

series

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Input power:	36W
Rated voltage:	100-240V
Rated frequency range:	50Hz / 60Hz
Storage temperature:	from 0oC to 45oC
Working temperature:	from 5oC to 35oC
Humidity during operation:	from 25% to 75%

#### Security:

Before using the product, read this manual and follow the recommendations contained therein when using the device.

#### General description [page 2- f1]:

1. Switch.
  2. Intensity indicator / LCD display.
  3. Replaceable lamp.
  4. Flash button.
  5. DC input.
  6. Lamp type (HR/hair removal, SR/skin rejuvenation, AC/acne treatment).
  7. Remaining impulse.
  8. Device temperature (light OK, flashing means the temperature is too high).
  9. Energy level.
  10. Fan status (light is OK, flashing means the temperature is too high).
- Additionally: storage case, charger.

#### Features [page 2- f1]:

Shine has three functions and separate lamps for each of them.

- Function 1: IPL hair removal  
 Function 2: IPL skin rejuvenation  
 Function 3: IPL acne removal
- A. Hair removal
  - B. Wrinkle reduction
  - C. Acne treatment

#### Contraindications:

1. The device is not intended for children under 15 years of age. Young people aged 15–18 may use the device with the consent and/or supervision of their parents or guardians. Adults over 18 years of age can use the device freely.
2. Do not use if the skin in the area to be treated is tanned or red. Using the device on tanned or red skin may cause temporary or even permanent discoloration of the skin.
3. Do not use if you are pregnant, trying to get pregnant, think you may be pregnant or are breastfeeding. The Shine device has not been tested for this purpose and therefore should not be used in these cases. Hormonal changes that occur during pregnancy or breastfeeding may increase the risk of skin damage or skin sensitivity.
4. Do not use the device if you suffer from diseases such as epilepsy, epileptic seizures or migraine. Energy flashes can cause severe headaches and migraines.
5. Do not use on areas of skin covered with tattoos, dark brown or black spots (such as dark freckles, moles, or warts). Tattoos or dark spots can absorb too much light energy, which can cause skin swelling, temporary or permanent blisters, discoloration, or skin whitening.
6. Do not use the device if you have a tendency to develop keloid scars.
7. Do not use on the scalp, ears, neck, nipples or mucous membranes (labia minora, vaginal and anal areas). In those parts of the bikini area where the skin is more pigmented, a lower light intensity level should be used. Use in these places may cause skin swelling, blistering or temporary discoloration or whitening of the skin.
8. Do not use if the treatment area has any active skin inflammation, including acne, psoriasis, eczema, burns, herpes, abrasions or rash. These factors can make the skin more sensitive to light and make the treatment uncomfortable or cause swelling, blisters or infections on the skin. Once your skin has healed sufficiently, consult your doctor before using the device.
9. Do not use the device if you are taking medications that make your skin sensitive to light, including non-steroidal anti-inflammatory agents (e.g. aspirin, ibuprofen, cetaminophen), tetracyclines, phenothiazines, thiazide, diuretics, sulfonyleureas, sulfonamides, DTIC, fluorouracil, vinblastine, griseofulvin, Alpha Hydroxy Acids (AHAS), Beta-Hydroxy Acids (BHAS), Retin-A, Accutane and/or topical retinoids.
10. Do not use if you have a disease associated with photosensitivity, e.g. porphyria, polymorphic photoeruption, solar urticaria, lupus, etc. These conditions make the skin more sensitive to light and the procedure may be unpleasant or painful. skin swelling, blisters or infections.

#### Handling Warnings:

1. The product will not turn on without inserting the lamp, please insert the lamp before connecting to the power. When the device is operating, do not cover the inlet and ventilation air openings, so as not to adversely affect heat generation.
2. Wear safety glasses when operating the device.
3. Before use, clean the area of skin to be treated. It is not allowed to apply any cosmetics to the skin immediately before the treatment. This may damage the device and affect the effect of the treatment.
4. Adjust your energy level from low to high. The higher the energy level, the better the result. Immediately after the treatment, skin redness may occur. Ice packs can be used for this purpose.
5. After phototherapy treatment, use a moisturizing preparation and preparations with high sun protection.

Table comparing skin and hair color:

✓ applies    ✗ not applicable

Skin color \ Hair color	White	Ivory shade	Beige	Light brown	Brown	Dark brown and black
	Black	✓	✓	✓	✓	✓
A dark brown	✓	✓	✓	✓	✓	✗
Brown	✓	✓	✓	✓	✓	✗
Golden	✓	✓	✓	✓	✓	✗
Red	✗	✗	✗	✗	✗	✗
White	✗	✗	✗	✗	✗	✗
Gray	✗	✗	✗	✗	✗	✗

More information about the company and products [page 2 -f2]:

- [www.beautify.eu](http://www.beautify.eu)
- [www.youtube.com/beautifyeu](https://www.youtube.com/beautifyeu)

6. After the treatment, use a damp cotton swab or cloth to clean the lamp from any remaining hair, sebum or sebum.
7. It is recommended that the device and accessories be kept away from water and moisture. Store the device in a ventilated and dry environment and keep it away from children.
8. It is prohibited to place metal objects in the heat emission hole. Otherwise there is a risk of electric shock.

Before starting the procedure:

1. Clean the cassette light output window with a damp cloth or cotton swab. This is done to remove hair and dust from the treatment window.
2. Hold the lamp above and below, then insert or pull parallel.
3. Connect the power cord to the Shine device.
4. Place the Shine device near an electrical outlet. Away from all kinds of liquids.
5. Make sure the lamp is installed correctly before turning on.
6. Once the power cord is plugged into the wall outlet, the status lights will slowly flash green.
7. Press the power button for more than 2 seconds to turn on/off the device's power. Then the fan starts and the LED lights up. Lamp type, remaining pulse, energy level (default level 1), thermometer icon, fan icon will be displayed on the LED [page 2- f1].
8. The status indicator flashes to indicate that the device is ready. When the work is finished, the indicator will turn off. When the treatment window touches the skin, the indicator turns green, which means the device is ready to emit light.
9. Press the power button for a while to adjust the energy level from low to high, level 1 is the lowest and level 5 is the highest.
10. The treatment window should directly touch the area of the body being treated.
11. The device will turn off automatically if not used within 5 minutes.

The device has two operating modes: automatic and manual. Manual mode turns on as the default mode. To start automatic mode, hold the flash button for a few seconds. If you want to return to manual mode, press the button again for a few seconds.

Perform a skin sensitivity test:

If you are using the device for the first time or on a new area of the body, select a small area near the area you intend to treat to perform a skin sensitivity test. Use the lowest setting of the device. Follow the treatment instructions below, but emit no more than 2-3 pulses. After the skin test, wait 60 minutes and check if your skin does not react. If after 60 minutes the skin appears normal and there is no redness or other skin reaction, you are ready for the full treatment.

Choose the right intensity:

Choose the appropriate intensity depending on your skin tolerance and increase it gradually from one to five levels. If the intensity is too weak and there are no skin reactions, you can gradually increase the intensity level.

Do not flash pulses in one place on the skin multiple times. This will create too much heat on the skin, which could cause skin swelling, blisters, or temporary or permanent changes in the color of your skin. After each flash, slowly move the device to an adjacent spot on the skin. Wait at least 20 seconds before returning to the previous area of skin.

Operation procedure:

Hair removal (depilation lamp)

Application:

- Mustache, armpit hair, body and leg hair, hairline hair on the forehead and bikini area.
- The effect varies depending on skin color. The device is not suitable for white and light colors

How to perform a hair removal treatment:

Usually it is about 7-9 treatments. Perform the first 3 treatments once a week. For the next 4-9 treatments, perform one therapy every 2-3 weeks. To maintain the effect of hair removal and to prevent the hair follicles from reactivating, repeat the treatments regularly to ensure that the skin remains smooth. If treatments are not continued, the hair follicles will become active again and the hair will begin to grow back. To maintain the results, simply repeat the treatments every 4-8 weeks to ensure that the hair follicles do not become activated and that hair regrowth remains inhibited. The intervals between treatments may vary depending on the individual rate of hair regrowth.

1. Before using the Shine device, prepare your skin by shaving the selected area to be treated. The skin should also be washed.
2. Wear safety glasses to avoid IPL light stimulating your eyes.
3. Select the appropriate intensity (see "select intensity") and start at a low level.
4. Allow the device window to touch the selected areas vertically and treat the areas one by one. Don't repeat the same place over and over or skip any area.
5. Start on the initial area and repeat step 4, 2-3 times.
6. After treatment, it is recommended to use moisturizing preparation after depilation.
7. Turn off the device and use a damp cloth or cotton swab to clean the treatment window.
8. Place your device and accessories in a case or box and store them in a safe place.

Skin photorejuvenation:  
(wrinkle reduction lamp)

The aim of rejuvenating selected areas is to use short and safe light pulses directed towards the skin surface. The light wave stimulates skin cells to produce collagen and elastin, which accelerate skin

regeneration. The therapeutic effects of IPL radiation include improved blood supply and oxygenation of the skin, an increase in the number of collagen and elastic fibers, better moisture absorption, and regulation of the pigmentation and keratinization process.

Effects:

- elimination of symptoms of skin aging
- reconstruction of collagen fibers
- firm and elastic skin
- shrinking of pores
- smoothing wrinkles
- improved skin hydration
- lightening discolorations

How often to perform a photorejuvenation treatment:

To obtain the desired effects, a series of treatments must be performed. It is recommended to carry out a series of 8 treatments, one treatment each week.

1. Before starting the treatment, you should thoroughly cleanse your skin of makeup and any cosmetics. You should also wipe the lamp.
2. Wear safety glasses to avoid IPL light stimulating your eyes.
3. Select the appropriate intensity (see "select intensity") and start at a low level.
4. Let the device window touch the selected areas vertically and treat the areas one by one. Don't repeat the same place over and over or skip any area.
5. Start on the initial area and repeat step 4, 2-3 times. After the treatment, it is recommended to use a strongly moisturizing preparation.
6. Turn off the device and use a damp cloth or cotton swab to clean the treatment window.
7. Place your device and accessories in a case or box and store them in a safe place.

Acne treatment:

(acne treatment lamp)

Acne treatment using IPL is intended for people with medium severity of acne lesions. IPL light destroys rapidly multiplying bacteria responsible for inflammatory changes in acne. The action is similar to antibiotics that are used topically or systemically. However, they can cause a number of side effects, damaging the liver and stomach.

Shine laser therapy works directly on the causes of acne. It acts on the so-called molecules synthesized by anaerobic bacteria. porphyrins. These make the bacterium particularly sensitive to light. When porphyrin is exposed to visible light, it becomes chemically active and reaches an excited state. This causes the formation of singlet oxygen, which combines with the cell membrane and destroys anaerobic bacteria responsible for the formation of acne lesions and sebaceous gland cells.

IPL light therapy must not be combined with antibiotic therapy. You should also wait a few weeks after completing such treatment.

How often to perform acne treatment:

The treatment period lasts approximately 8 - 16 weeks (2-4 months), and treatments are performed every 1-2 weeks. It is best to perform two irradiations of all affected areas during the procedure, approximately 10 minutes apart. Starting from the second week, the first signs of improvement are visible (the degree of inflammation in the affected areas and a slight reduction in the symptoms themselves occur).

After the entire series of treatments:

Acne spots are almost completely eliminated. However, those that remain - no longer active - disappear over the next weeks. After healing, it is recommended to use skin photorejuvenation therapy, which is able to stimulate collagen production in order to reduce the number of acne scars. During the entire treatment, it is recommended to use appropriate care cosmetics intended for oily and acne-prone skin.

How to perform an anti-acne treatment:

1. Before starting the treatment, you should thoroughly cleanse your skin of makeup and any cosmetics. You should also wipe the lamp.
2. If inflammation occurs in the area to be treated, the skin should be healed before using the device.
3. Wear safety glasses to avoid IPL light stimulating your eyes.
4. Select the appropriate intensity (see "select intensity") and start at a low level.
5. Let the device window touch the selected areas vertically and treat the areas one by one. The interval between pulses should be at least 20 seconds.
6. After the treatment, it is recommended to wash your face and use a moisturizing preparation.
7. Turn off the device and use a damp cloth or cotton swab to clean the treatment window.
8. Place your device and accessories in a case or box and store them in a safe place.

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**[EN] Rules for handling electrical and electronic equipment:** Electrical and electronic equipment as well as batteries and accumulators must not be combined with other household waste. This is indicated by the crossed-out wheeled bin symbol on the packaging. Dispose of waste electrical and electronic equipment as well as batteries and accumulators at your local separate collection point or to your dealer. Electrical and electronic equipment and batteries and accumulators may contain hazardous substances, mixtures and components that are harmful to the environment and human health. Incorrect handling or damage to them may lead to damage to health or environmental pollution when the equipment is further disposed of or recycled. The device was placed on the market after August 13, 2005

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**[EN] Summary declaration of conformity:** The manufacturer hereby declares under his sole responsibility that the above product complies with Directive 2014/30/EU of the European Parliament and of the Council of 26 February 2014 on the harmonization of legislation. Member States in the field of electromagnetic compatibility and Directive 2011/65/EC of the European Parliament and of the Council of 8 June 2011 on the restriction of the use of certain hazardous substances in electrical and electronic equipment and Commission Delegated Directive (EU) 2015/863 of 2011. On 31 March 2015, amends Annex II to Directive 2011/65/EU of the European Parliament and of the Council on the list of restricted substances.

Model: Shine Blush